

## Customer and Patient Details

Patient Name/ID/Ref: \_\_\_\_\_

Request for Quote  Purchase Order

Contact Person: \_\_\_\_\_

PO Number: \_\_\_\_\_

Email: \_\_\_\_\_

Phone: \_\_\_\_\_

Ship to: \_\_\_\_\_

Bill to: \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

Post Code: \_\_\_\_\_

Post Code: \_\_\_\_\_

## Basic Seat Frame Size

(Includes: Frame, Arm Rests, Seat Cushion, Rear Feet)

Goodwood 2	Goodwood 3	Goodwood 4	Goodwood 5	Goodwood 6
<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Seat Width: 27 cm Seat Depth: 31 cm Seat Height: 25 cm	Seat Width: 33 cm Seat Depth: 35 cm Seat Height: 31 cm	Seat Width: 38 cm Seat Depth: 40 cm Seat Height: 38 cm	Seat Width: 43 cm Seat Depth: 46 cm Seat Height: 45 cm	Seat Width: 48 cm Seat Depth: 50 cm Seat Height: 51 cm
Arm Rest Height: <input type="text"/> cm Measured from the floor to top edge of the arm rest				

## Additional Padding Options

With Reduction Pads  ½"  1"  2"  3"  Tapered

---

Back Pad  ½"  1"  2"  3"  With sacral lump  
Not available with 3" pad

---

Contoured Back Pad  1"  2"  3" Number of 15mm spacers to move back forward on chair (1- 4)

---

Comfort Pad for Footplate  ½"  1"  2"  3"

## Upper Body Positioning Options

Corner Back:  Number of 15mm spacers to move back forward on chair (1- 4)

PivotFit Harness	Stayflex	Monoflex Underarm Release
<input type="checkbox"/> Standard <input type="checkbox"/> or Dynamic	<input type="checkbox"/> With Zip <input type="checkbox"/> Without Zip	<input type="checkbox"/> XS <input type="checkbox"/> Sm
<input type="checkbox"/> XXS <input type="checkbox"/> XS <input type="checkbox"/> Sm	<input type="checkbox"/> Sm <input type="checkbox"/> XS <input type="checkbox"/> Sm	<input type="checkbox"/> XS <input type="checkbox"/> Sm

## Pelvic Positioning and Support

### Pelvic Positioning Belts

Pad Size:  Extra Small (32 cm pad width)  Small (38 cm pad width)

Buckle Type:  Push Button  With Security Cover  
 Swivel

### Knee Abduction Pad (Pommel):

6 cm diameter  8.5 cm diameter

Note: Comes with Quick Release Mount

Groin Straps  Velcro Strap

## Lower Body Positioning

Foot Straps:  48 cm Long

Foot Shoes:  Small  Large

Ankle Huggers:  XS  Small

Velcro  Side Release

## Mobile Options

Mobile Base  Skis with Casters

Millie Mova  Mobile Frame

Push Handles  Colibri Hi Lo Frame  
 (Only for Goodwood 2 & 3)

## Additional Optional Accessories

Arm Rest Only  Heel Bar  Foot Board  Stabilising Skis

Cross Grab Rail:  Full Length  ¼ Length

Head Rest:  Flat Only  With Wings

Tray  With Tray Securing Straps

## Upholstery Colour

Standard Material:  Red  Blue  Time Square Taupe Faux Leather

### Vinyl

Apple Mint <input type="checkbox"/>	Cobalt Black <input type="checkbox"/>	Lupin Blue <input type="checkbox"/>
Cool Blue <input type="checkbox"/>	Eden Green <input type="checkbox"/>	Marigold <input type="checkbox"/>
Candy Pink <input type="checkbox"/>	Grape <input type="checkbox"/>	Pillar Box Red <input type="checkbox"/>
Blossom Pink <input type="checkbox"/>	Jasmine White <input type="checkbox"/>	Deep Rose <input type="checkbox"/>
Cocoa <input type="checkbox"/>	Lilac <input type="checkbox"/>	Ultra Violet <input type="checkbox"/>
Rouge <input type="checkbox"/>	Sapphire <input type="checkbox"/>	

### Faux Leather

Almond   
 Greenwich Green

Extra charge for non - standard colours

## Additional Notes

## Dealer Stamp



DCS Joncare Limited  
 7 - 8 Radley Place • Radley Road Industrial Estate • ABINGDON • Oxfordshire • OX14 3RY  
 Tel: 01235 523353 Fax: 01235 531 019 info@joncare.co.uk www.joncare.co.uk  
 Company Registration No. 2647496 VAT Registration No. 559160526 © Joncare 2016 EOE



JLFT04162