

The **Stealth** Tablet Mount

Tablets are used by wheelchair users as communication and learning tools, as well as, to stay connected. To be able to use these devices it is important to have a mounting system that makes them accessible and reachable.

The **TWBM480-TAB** fits virtually any tablet (compatibility includes but is not limited to iPad 1, iPad 2, iPad 3, iPad4, HP TouchPad, Asus EEE Pad Memo, Barnes & Nobles NOOKcolor, BlackBerry PlayBook, Dell Streak 7, HTC Flyer, Samsung Galaxy tablets and more). It includes adjustable tablet mounting system for wheelchair - TWB480-17FMD with **Arms240 Mount** and **2 Ratchet Handles**. To make adjustment easy, Stealth has swapped normal screws on the links of the system with ratchet handles.

This mounting system can be mounted under armrests, and is available in a swing-away option, so that it can be swung away for easy patient transfers, and clicked back to its functional position.

