

Millie-Mova

“The single best invention for safe manual handling of people since the hoist”

Gavin Wright, Moving & Handling Consultant

Scan the QR Code to see a movie of the Millie-Mova in action



An easy to use, safe manual handling tool - the **Millie-Mova**

The Millie-Mova is a unique, easy to use, and safe manual handling system designed to help a carer move an elderly or disabled person safely while seated in a chair. The Millie-Mova attaches to the base of any four-legged chair, and has a spring mechanism that allows the carer to engage the wheels to glide a chair from one place to another with an occupant in place, without risking possible injuries to their back, shoulder, knees, or hips while moving the person.

Unlike other solutions to the problem of moving elderly or disabled people whilst safely seated, a chair fitted with the Millie-Mova retains high stability, and incorporates an automatic brake, when unoccupied. This means that a person unsteady on their feet can lean on the chair without it sliding away from them, which otherwise could lead to a fall.

► Millie-Mova installed on a chair (chair not included)



Benefits of the **Millie-Mova**



Easy to use



Simple to install



Fits most 4-legged chairs



Prolongs chair and floor life



Improves dignity and safety of the client in the chair



Prevents carer from risk of injury caused to the back, hip, knees and shoulder while moving a client's chair

Millie - Mova

| Part Number | Description |
|----------------------|---|
| Millie-Mova - 40 | For chairs up to 40 cm wide |
| Millie-Mova - 47 | For chairs up to 47 cm wide |
| Millie-Mova - 40M | For chairs with metal legs up to 40 cm wide |
| Millie-Mova - 47M | For chairs with metal legs up to 47 cm wide |
| Millie-Mova - 40PED | For chairs up to 40 cm wide. Includes a pedal for easier operation of the Millie-Mova. |
| Millie-Mova - 47PED | For chairs up to 47 cm wide. Includes a pedal for easier operation of the Millie-Mova. |
| Millie-Mova - 40PEDM | For chairs with metal legs up to 40 cm wide. Includes a pedal for easier operation of the Millie-Mova. |
| Millie-Mova - 47PEDM | For chairs with metal legs up to 47 cm width. Includes a pedal for easier operation of the Millie-Mova. |