



Supporting your needs

# The Powerchair Assessment Game Changer



Assessments made simple with  
the i-Drive kit and Loonz game

## The Powerchair Assessment Game Changer



### The i-Drive assessment kit combined with the Loonz game has revolutionised the powerchair assessment and prescription process, ensuring better outcomes.

The i-Drive kit gives you the ability to **assess individuals multiple times in their existing seating system without the need for a powerchair.**

This means you can gain a **true assessment of an individual's functional ability** as they can focus on operating the switches rather than struggle with an unfamiliar seating system.

The i-Drive in combination with the Loonz game ensures individuals are comfortable with power controls, helping to build confidence within the

user. It can also help build confidence in parents who may be nervous about their child operating a heavy powerchair.

With the Loonz game the individual is assessed via a series of interactive levels. Here, they control Tiki the elephant explorer, helping him collect coins while avoiding obstacles and dangers as he flies in his hot air balloon. As the levels progress, the difficulty increases, **helping individuals improve their hazard awareness.**

With the Loonz game you can **record progress and make informed decisions before issuing a suitable powerchair and power controls.** It can assess for a range of controls from head arrays, to sip and puff, tray controls, and toe controls.

The i-Drive kit and the Loonz game are ground-breaking pieces of technology that enable therapists to carry out true powerchair assessments, while adding the element of fun for the user.



## Select your i-Drive Assessment Kit

### Professional Kit

For those who require various drive controls to assess for a whole range of conditions. The Professional Kit enables healthcare professionals to make complex decisions, and comes with a comprehensive selection of switches to cover all switch access points (head, hands, etc.).

### Head Array Kit

For those assessing people with a degree of head control and cognitive ability. The Head Array Kit is the perfect assessment kit for assessing those where movement is more prominent in the head than any other part of the body.

### Switch Kit

For those assessing first-time powerchair drivers. The Egg Switch components can be placed at various switch access points (head, hands, etc.) to assess the most appropriate placement for the individual.

## Select your Stealth I-Drive Kit

Code	Description	Professional	Head Array	Switch
IDHBT300-1	i-Drive Head Array with Bluetooth	✓	✓	
IDTF100-2	Eclipse Tray	✓		
IDM-MICRO-9	Mo-Vis Micro Joystick	✓		
IDM-MULTI-9	Mo-Vis Multi Joystick	✓		
IDPMP-1	Precision Mini Portable Joystick (PMPJ)	✓		
IDM-ALLROUND-9	Mo-Vis All Round Joystick	✓	✓	
IDHBT501	Bluetooth Proximity Sensor	✓		✓
IC12A-BC	12V Battery & Charger for i-Connect	✓	✓	✓
P002-31	Joystick Interface - Omni DX	✓	✓	
ESRED	Egg Switch - Red		✓	✓
ESGRN	Egg Switch - Green		✓	✓
ESM	Egg Switch Mount		✓	
SU264	13" Lateral Rod		✓	
IDHBT500	i-Drive Smart Box with Bluetooth			✓
SK-BOX	Stealth Kit Box - Small	✓	✓	✓

# i-Drive Professional Kit

## i-Drive Kit Components



Code: IDHBT300-1  
i-Drive 4.0 Head Array with Bluetooth



Code: IDHBT501  
i-Drive Proximity Sensor



Code: IDTF100-2  
Eclipse Tray



Code: IC12A-BC  
12V Battery & Charger for i-Connect



Code: IDM-ALLROUND-9  
Mo-Vis All Round Joystick



Code: IDM-MICRO-9  
Mo-Vis Micro Joystick



Code: IDM-MULTI-9  
Mo-Vis Multi Joystick



Code: IDPMP-1  
PMPJ Joystick



Code: P002-31  
Joystick Interface - Omni DX



Code: ESRED  
Egg Switch - Red



Code: ESGRN  
Egg Switch - Green



Code: ESM  
Egg Switch Mount



Code: SU264  
13" Lateral Rod



Code: IDHBT500  
i-Drive Smart Box with Bluetooth

Small Stealth Kit Box included FOC with the Professional Kit, and included within the price of Head Array and Switch Kit.