


Pelvic Positioning

Standard to describe Pelvic Positioning Belts



COMPATIBLE WITH
BS 8625

Webbing Width	Pad Length	Pad Width
S = 1" (25 mm)	32 cm	4 cm
M = 1.5" (38 mm)	36 cm	5 cm
L = 2" (51 mm)	38 cm	6 cm
	46 cm	7 cm
	62 cm	9 cm
	82 cm	10 cm
	92 cm	



Pelvic Positioning Belts

Using a flexible measuring tape, measure the distance over the body - this corresponds to the belts padded length.

Where to fix a belt to correct for posterior tilt



ISO 16840-3:2022 Wheelchair seating — Part 3: Determination of static, impact, and repetitive load strengths for postural support devices
BS 8625:2019 Selection, placement and fixation of flexible postural support devices in seating

Hip Belts 4-Point

	Webbing Width	Pad Length	Metal Push Button	Standard Push Button	Rehab Latch	Plastic Side Release	Swivel Buckle
Centre Pull	S (1")	32 cm	132453*			101749	
	S (1")	38 cm	109936*			109448	108090
	M (1.5")	36 cm		100131	101527	101121	
	M (1.5")	46 cm		109240*	109851	109400	107987
Dual Pull	L (2")	62 cm		109202**	109783	109318	
	S (1")	32 cm				101800	
	S (1")	38 cm				109622	
	M (1.5")	36 cm				101459	
Rear Pull	M (1.5")	46 cm				109585	
	L (2")	62 cm				109547	
	S (1")	38 cm	110055*			109707	
	M (1.5")	46 cm		111397*	111670	111519	
	M (1.5")	62 cm		109271*	109899	109660	

*Also accepts the 165635 reduced access or 100353 security covers. **Also accepts 175849 reduced access or 100599 security covers. *Small Metal Push Button, accepts 176525 Security Cover.

Leg Harness

Size	Code	Measuring for a Leg Harness
Small	141325	<p>Measure from the Anterior Superior Iliac Spine (ASIS) to the point where the inner thigh touches the seat.</p> <p>S 8"-11" (20-28 cm) M 11"-14" (28-36 cm) L 14"-16" (36-41 cm)</p>
Medium	141561	
Large	141707	

Hip Belts 2-Point

	Webbing Width	Pad Length	Metal Push Button	Standard Push Button	Rehab Latch	Plastic Side Release	Swivel Buckle
Centre Pull	S (1")	32 cm	132538*			100087	100254
	S (1")	38 cm	109066*			108557	199753
	M (1.5")	36 cm		100070	101336	100100	101176
	M (1.5")	46 cm		108281*	108946	108489	199814
Plus Size Centre Pull	L (2")	62 cm		108243**	108908	108403	
	L (2")	82 cm		199876**			
Dual Pull	L (2")	92 cm		199937**			
	S (1")	32 cm				101695	
	S (1")	38 cm				108717	
	M (1.5")	36 cm				101282	
Rear Pull	M (1.5")	46 cm				108632	
	L (2")	62 cm				108595	
	S (1")	38 cm	109165*			108861	
	M (1.5")	40 cm		110987*	111205	111090	
	M (1.5")	46 cm		108328**	109035	108793	
Evoflex®							
Centre Pull	S (1")	38 cm	105211*				100179
	M (1.5")	46 cm		105013*	191061		100513
	L (2")	62 cm		104979**	191443		

*Also accepts the 165635 reduced access or 100353 security covers. **Also accepts 175849 reduced access or 100599 security covers. *Small Metal Push Button, accepts 176525 Security Cover.

Hook and Loop Belts

	Webbing Width	Pad Length	Metal Push Button	Standard Push Button	Rehab Latch	Plastic Side Release
Hook and Loop Belts 2pt						
Centre Pull	S (1")	38 cm	110444*			
	M (1.5")	46 cm		110130*	110369	
	L (2")	62 cm		110093**	110291	
Rear Pull	M (1.5")	46 cm		110222*	110406	
Hook and Loop Belts 4pt						
Centre Pull	S (1")	38 cm	110901*			
	M (1.5")	46 cm		110604*	110826	
	L (2")	62 cm		110529**	110789	
Dual Pull	S (1")	38 cm				110741
Rear Pull	M (1.5")	46 cm		110659*	110864	
Sub ASIS Pads (for Hook and Loop Belts)						
Small	149031			Sub ASIS pads greatly enhance stability and improve control of obliquity and rotation while minimising abdominal compression.		
Medium	148973					

*Also accepts the 165635 reduced access or 100353 security covers. **Also accepts 175849 reduced access or 100599 security covers. *Small Metal Push Button, accepts 176525 Security Cover.

Hardware Options

Framesaver™ Clamp		Seat Tube Mounting Bracket Kit	
Frame Clamp		Cinch-Mount™ End Fitting	
Flat-Mount™ & Steel Sliders		Band Clamp	
Belt Extension Kit		Hardware Kit (191207)	

Reduced Access Cover

Size	Code	
Medium Belts	165635	(14 mm hole)
Large Belts	175849	
Security Cover		
Small Metal PB	176525	
Medium Belts	100353	
Large Belts	100599	

Buckle Style




Pull Style




Upper Body Positioning

Standard Stayflex				Monoflex				
	Size	With Zip	Exc Zip		Size	Underarm Release with Flat Mount	Centre Release with Flat Mount	Centre Release
	Small	103149	103248		Extra Small	114039	114107	114183
	Medium	103118	103224		Small	113858	114220	114268
	Large	103088	103194		Medium	113735	114305	114343
	Extra Large	103163	103279+		Large	113650	114428	114466
				Extra Large	113957	114565	114602	
Narrow Stayflex				Trimline Shoulder Harness				
	Size	With Zip	Exc Zip		Size	Standard		
	Medium	103354	103422		Small	102746		
	Large	103330	103392		Medium	102722		
	Extra Large	103378	103453		Large	102692		
PivotFit Shoulder Harness				Harness Hardware				
	Size	Standard	Dynamic		Plastic Flat-Mount™ & Sliders	Mounting Adapter	Strap Guides	Plastic Cinch Mounts
	XXS Small	102906	103064					
	XS Small	102876	103040					
	Small	102821	102999					
	Medium	102807	102968					
	Large	102777	102937					
Extra Large	102852	103026						

How to measure for a Shoulder Harness



How to measure for a Monoflex and Chest Strap



Lower Body Positioning

Ankle Huggers			
Size	Hook and Loop/ ErgoLatch™		How to measure for Ankle Huggers
Extra Small	102401		Measure circumference above ankle
Small	102371		XS 5.5"-7" (14-17 cm)
Medium	102364		S 6.5"-8" (17-20 cm)
Large	102357		M 7.5"-9" (19-23 cm)
Extra Large	102395		L 8.5"-11.5" (22-29 cm)
			XL 11"-13" (28-33 cm)
Calf Panel, Calf Strap, Toe Strap			
Size	Calf Panel	Calf Strap	Measuring for Calf Panels and Calf Straps
Small	169213	164027	For proper chair fit, measure chair width between tubing.
Medium	169145	163723	S 12"-15" (31-38 cm)
Large	169077	163419	M 15"-19" (38-48 cm)
Toe Straps	184681		L 19"-22" (48-56 cm)

Go Anywhere Positioning

Universal Elastic Strap				
	Size	Strap	Strap Extension	Length: 56" (Extensions available)
	Small	180553	182397	Width:
	Large	180799	182465	S 3" (8 cm)
	Extra Large	101411	181789	L 4.5" (11 cm)
				XL 6" (15 cm)
Rapid Dry Belts				
	Type	Size	Order Code	Scan for measuring info:
	One Piece with Grommets	Large	183073	
		Extra Large	182915	
	One Piece with Plastic Slides	Large	184209	
		Extra Large	184131	
	Two Piece with Grommets	Large	159795	
		Extra Large	165321	